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GB – ANTI-SNORING BANDAGE – 008227

The anti-snoring bandage keeps the mouth closed and prevents the lower jaw from moving backward during sleep. It can help reduce snoring by promoting nasal breathing.

Instructions for use: Place the bandage over your head and under your chin (see illustration) and thread your ears through the openings. Adjust the circumference using the adjustable Velcro fasteners so that the bandage fits comfortably and securely without causing pressure points. Do not fasten the bandage too tightly to avoid restricting blood circulation.

Safety instructions: If you suffer from sleep apnea or breathing difficulties, consult a doctor before using the bandage. Pregnant women and individuals with medical conditions should also consult their doctor before use. Do not use the bandage if your nose is congested or if you have respiratory conditions that make it difficult to breathe through your nose. If skin irritation or allergies occur, stop using the bandage and consult a doctor. This product is not suitable for children. Keep it out of the reach of children.

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