DE – Gebissreiniger

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EN - Denture cleaner

Mount two 1,5V round batteries into the battery compartment on the bottom side of the unit. Open the upper flap lid and put the denture into the (removable) clear plastic tray. Fill the tray with clear, fresh water and add your favourite denture cleaning agent. Shut the flap lid and activate the device by sliding the switch on the front into the "ON" position. Now the constant vibration, in conjunction with the cleaning agent, will clean your denture properly and thoroughly. After 6 - 10 minutes, you can switch off the device and take out your squeaky clean denture. Rinse the denture with warm water after each cleaning procedure. Cleaning: Wipe the Denture Cleaner with a moist cloth and wipe it dry afterwards. Don't dip the device into water.

Disposal of batteries: Don't put batteries into normal household litter. As a consumer, you are legally bound to return used batteries accordingly. You can use either public collecting points or any shop that sells the specific kind of battery that you want to return.

Disposal of used product: At the end of its lifespan, don't dispose this product/unit together with normal household litter. Bring unit to a public collecting point for the recycling of used electric and electronic devices. The referring symbol can be found either on the product, the instruction manual or the packaging.

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GHZ Matra AG, Max-Volmer-Str.14, 40724 Hilden, Germany