# HydroBlast PRO™

## Kabellose Munddusche

DE .....	4
EN .....	8
IT .....	12
FR .....	16
NL .....	20
HU .....	24
PL .....	28
RO .....	32

Abbildung 1 | Figure 1 | Figura 1 | Image 1 | Afbeelding 1 | Ábra 1 | Zdjęcie 1 | Ilustrația 1

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## **EN** WARNING

- Read the instruction before use. Keep the instruction in a safe place in case you need them again.
- When the HydroBlast Pro is not in use keep it dry and store it away from fluids.
- Do not use HydroBlast Pro in a shower or bath.
- Do not immerse the HydroBlast Pro in water.
- Do not use the HydroBlast Pro if either the appliance or accompanying accessories are damaged.
- Do not use if the HydroBlast Pro has been dropped or is not operating correctly.
- Do not aim the water stream under your tongue, in your ear or at your eyes.
- For use on gums only. The pressure of the HydroBlast Pro can damage delicate skin.
- Please consult your doctor or dentist before using the HydroBlast Pro.
- Only use recommended accessories with the HydroBlast Pro.
- Do not insert any foreign object into the HydroBlast Pro
- Keep the HydroBlast Pro and all its accessories away from heat.
- Do not use the HydroBlast Pro outside nor in the vicinity of aerosols or where oxygen is being administered.
- Do not under any circumstances use iodine, bleach or water insoluble concentrated oils with or inside the HydroBlast Pro.
- Do not fill the tank with anything other than water or your preferred mouthwash, unless otherwise directed by your doctor or dentist.
- Do not use HydroBlast Pro whilst wearing jewellery inside your mouth.
- If you are taking antibiotics or any oral hygiene medication, consult your doctor before using

the HydroBlast Pro.

- Do not use the HydroBlast if you have an open wound or other irritation inside your mouth.
- Children should **ONLY** use the HydroBlast Pro under strict adult supervision.

## PRODUCT DESCRIPTION

**HydroBlast Pro** is a portable, battery-operated cordless mouth water flosser. The HydroBlast Pro releases a fine stream of water for dental cleaning. The appliance comes with an interchangeable brush and water stream tips. The large volume of water used allows easy and effective dental cleaning at home or when travelling. It is safe for use on crowns, braces and other dental aids. Fill with water or your preferred mouth wash and choose the low or high pump setting.

**Intended use:** for oral rinsing and dental cleaning. Particularly effective when used on braces.

## HydroBlast Pro components

(See Fig. 1):

- A** Standard water pressure tip (blue mark)
- B** Refill opening
- C** On/off switch and speed control
- D** Easy-Grip handle
- E** Removable reservoir
- F** Brush attachment
- G** High water pressure (red mark)

## INSERTING BATTERIES

The battery compartment is located on the underside of the HydroBlast Pro. This can be opened and shut by turning using your fingernail or small coin, as illustrated. Requires 2 (not supplied) “AA”-batteries. (See Fig. 2)

## REMOVING AND REPLACING THE RESERVOIR

To remove the appliance’s reservoir, hold the

appliance in your hand and pull the plastic container downwards. Make sure that the silicone O-ring is attached to the reservoir. Align the container with the guide markers and slide upwards until it is seated firmly in place. **(See Fig. 3)**

## CLEANING & MAINTENANCE

Never immerse the appliance in water. Handwash only. Clean the appliance before using for the first time. Wash the container and tips in warm, soapy water and rinse with clean water. Always switch off the HydroBlast Pro after each use. Discard any remaining liquid. Clean the water container in warm, soapy water and rinse with clean water. The reservoir should be cleaned once a week. The tips should be cleaned regularly and replaced when worn.

## OPERATING INSTRUCTIONS

### To fill:

Ensure that the reservoir is attached properly. Open the top flap of the appliance. Fill with fresh water. Make sure that the silicone O-ring is fixed properly in place on the flap. Close the flap after filling **(See Fig. 4)**

**Advisory:** for even better results and ultimate freshness we recommend using mouthwash.

### Use:

Select the most suitable tip.

- Tip for standard water pressure (blue mark)
- Tip for high water pressure (red mark)
- Brush attachment

Insert the tip into the appliance. **(See Fig. 5)** Place the tip on your gum between two teeth. **Close your lips gently over the tip until it has reached the correct position to avoid any splashing.** Lean over a washbasin and set the appliance to the LO (low) setting. Allow the water to flow from the mouth into the washbasin.

Move the switch to position 2 to operate on HI (high). Move the tip along the gum until you can feel it between the next two teeth.

### Teeth cleaning tips

- The tip should be aimed at a 90-degree angle towards the gum line.
- Begin with the back teeth, working forwards. Aim the water stream between the teeth and gum line.

### Orthodontic cleaning

The HydroBlast Pro is effective at removing food particles, which are trapped in braces by using a combination of the water stream and brush.

## TROUBLE SHOOTING

### ▪ Water is leaking from the appliance?

Ensure that the silicone O-rings are affixed properly to both the reservoir opening and water inlet.

### ▪ No water is flowing through the tips?

- A. The tip may be clogged. Check and replace if necessary.
- B. The water container may be empty.
- C. The feed tube in the water container may be blocked.

### ▪ Pump pressure is too weak?

The batteries may need replacing.

### ▪ Water is spraying too heavily from your mouth?

When using the appliance, please ensure that you keep your lips closed around the tip as far as possible and select the lowest power setting at the start of use.

■ Please dispose of the appliance in an environmentally friendly manner. Do not dispose of the appliance with general household waste. Take it to a recycling centre that accepts electrical and electronic waste products. You can obtain further information from your local council.

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CH: MediaShop AG, Industriering 3, 9491 Ruggell, Lichtenstein
EU: MediaShop Holding GmbH, Schwarzottstraße 2a, 2620 Neunkirchen, Austria